Duraheal

- GB** Transparent Pericardium Membrane for Duraplasty
Instructions for Use
- DE** Transparente Perikard-Membran für die Duraplastik
Gebrauchsanweisung
- FR** Membrane péricardique transparente destinée à la plastie durable
Notice d'utilisation
- ES** Membrana transparente de pericardio para duraplastia
Información de uso
- IT** Matrice di collagene naturale per la rigenerazione della dura madre
Istruzioni per l'uso
- PT** Membrana pericárdica transparente para enxerto dural
Instruções de uso
- TR** Dura-plasti için saydam perikard membran
Kullanım talimatları
- EL** Διαφανής περικαρδιακή μεμβράνη για duraplasty
Οδηγίες χρήσης

Manufacturer

MBP Medical Biomaterial Products GmbH
Lederstraße 7 / 19306 Neustadt-Glewe
info@mbp-gmbh.de / www.mbp-gmbh.de

CE0483

Rev 09/2021-03

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DESCRIPTION

DURAHEAL is a watertight, foil-like collagen matrix manufactured from porcine pericardia (the pericardium is a membranous sac that surrounds the heart). The pericardia are sourced from establishments regularly inspected by veterinary authorities. In a validated procedure, the pericardia are freed from cells, cellular constituents and substances stored in the tissue such as fat and are then sterilised with ethylene oxide. This preservation ensures that the biomechanical properties of the original tissue are very largely retained and the surface is not cross-linked.

DURAHEAL is fully resorbable. DURAHEAL is available in a variety of sizes.

DURAHEAL is supplied dry and sterile in a double-peel package and contains no preservatives. The surgical membrane is a strong and biocompatible implant that is incorporated into the recipient tissue by cellular and microvascular infiltration.

COMPOSITION

Native, sterile, acellular collagen matrix made of porcine pericardium.

INDICATIONS

For use during neurosurgical operations as a dura substitute for the watertight repair of dura mater defects. DURAHEAL is indicated whenever a primary wound closure is impossible or appears unsuitable for other reasons.

CONTRAINDICATIONS

- DURAHEAL is only designed for the indicated use and is only sold for this purpose.
- DURAHEAL should not be used for patients with a known or suspected hypersensitivity to porcine collagen.

PREPARATION OF DURAHEAL BEFORE IMPLANTATION

- Open the cardboard box and remove the sterile double-peel package outside the sterile field and inspect for damage.
- Provided the packaging is intact, remove the inner sterile package from the outer one and deliver to the sterile field.
- After opening the inner sterile package, use sterile, atraumatic forceps to grasp the edge of DURAHEAL and remove it using the aseptic technique.

Transparent Pericardium Membrane for Duraplasty

- DURAHEAL must be rehydrated by soaking it in isotonic saline in a sterile dish for at least 5 minutes.

Clinical experience has shown that antibiotics capable of tissue penetration, such as gentamicin, may be added. The long-term effects of antibiotic treatments on DURAHEAL have not been investigated.

Antibiotics must only be used in accordance with the manufacturer's instructions.

- Keep the patch immersed in sterile saline until ready to use. **THE PATCH MUST BE KEPT MOIST AT ALL TIMES.**

IMPLANT INSTRUCTIONS

The patch can be trimmed during the operation to meet the surgeon's needs.

- Inspect both sides of the DURAHEAL patch. If one side appears smoother, implant the smoother side so that it faces the neural tissue. The limited thickness of the film-like collagen matrix enables the overlay technique to be used.
- DURAHEAL can be attached to the edge of the host tissue by sutures, clips, staples or fibrin glue.

- When implanting by suture, the suture bites must be placed 2-3 mm from the edge of the graft.

STORAGE

DURAHEAL is a sterile medical product that must be stored protected from direct heat and sunlight in a clean and dry place at room temperature.

Store in the original packaging.

The product must not be used after the expiry date.

WARNINGS

Do not resterilise.

Do not use if the packaging has been opened or damaged.

Do not use if the seal is broken or damaged.

After use, handle and dispose of all unused products and packaging materials in accordance with recognised medical practice and the applicable national and regional environmental protection regulations for the disposal of packaging and biological materials.

DURAHEAL must not be re-used once the packaging has been removed and/or

Transparent Pericardium Membrane for Duraplasty

has come in contact with a patient, since there is then an increased risk of contamination with an associated risk of infection.

USE IN PREGNANCY AND LACTATION

No studies have been conducted on the use of DURAHEAL in pregnancy or lactation or on its effects on human fertility. Therefore, before use the surgeon must carry out an individual assessment of the benefits for the mother and the possible risks to the child.

PRECAUTIONS

Avoid contact of DURAHEAL with powdered surgical gloves; if necessary these must be rinsed beforehand. Before surgery, patients or their representatives must be informed about the complications that could arise from the use of this product. As with any surgical procedure, these include the possibility of dehiscence at the surgical site and infection.

SIDE EFFECTS

No side effects have been reported to date. However, a reaction to xenogenous

collagen – although extremely rare – is possible. In the case of patients with an allergic predisposition or with a known intolerance to porcine collagen, emergency measures must always be readily to hand.

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