Xenoderm

GB Temporary biological skin replacement / Instructions for use
DE Temporärer biologischer Hautersatz / Gebrauchsanweisung
FR Remplacement cutané biologique temporaire / Notice d'utilisation
ES Sustituto cutáneo biológico temporal / Información de uso
IT Trapianto di cute naturale provvisorio / Istruzioni per l'uso
PT Substituição biológica temporária da pele / Instruções de uso
TR Geçici biyolojik deri ikamesi / Kullanım talimatı

CE0483

**Product overview / Produktübersicht / Aperçu du produit /
 Resumen del producto / Panoramica del prodotto /
 Visão geral do produto / Ürün özeti**

REF	Size/Größe/ Dimension/ Tamaño/ Dimensioni/ Tamanho/ Boyut	Content/Inhalt/ Contenu/ Contenido/ Contenuto/ Conteúdo/ İçerik	PZN	Product UDI-DI/Produkt UDI-DI/ UDI-DI du produit/ UDI-DI del producto/ UDI-DI del prodotto/ UDI-DI do produto/ Ürün UDI-DI
0505X	5*5 cm	10	2738218	04260230973017
1010X	10*10 cm	5	2738388	04260230973024
1020X	10*20 cm	5	2738394	04260230973031
1030X	10*30 cm	5	2738425	04260230973048
1040X	10*40 cm	5	–	04260230973055
0505XZ	5*5 cm	1	–	04260230973017
1010XZ	10*10 cm	1	–	04260230973024
1020XZ	10*20 cm	1	–	04260230973031
1030XZ	10*30 cm	1	–	04260230973048
1040XZ	10*40 cm	1	–	04260230973055

Serious incident reporting / Meldung schwerwiegender Vorkommnisse / Déclaration d'incident grave / Notificación de incidente grave / Segnalazione di incidenti gravi / Comunicação de incidente grave / Ciddi olay bildirimi

In the event of a serious incident involving Xenoderm, please report it to MBP Medical Biomaterial Products GmbH and to the competent authority in your Member State.

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URL: <http://www.mbp-gmbh.de/xenoderm>

English

Xenoderm – Temporary biological skin replacement

Definition

Xenoderm acts as a temporary skin replacement within the functional scope of the patient's own skin in order to:

- reduce fluid and heat loss
- promote the growth of epithelial cells and granulation tissue
- protect granulating wounds.

Composition

Xenoderm is a sterile, temporary biological collagen matrix of porcine origin.

The product consists of collagen obtained from porcine skin.

The main component relevant to the intended principal mode of action is collagen.

The collagen content is more than 80% based on dry weight.

No other components relevant to the intended principal mode of action are contained.

The biological safety of Xenoderm has been assessed as part of the biological evaluation for the intended application. During the course of wound healing, the collagen matrix is partially biologically resorbed in the wound area.

Indications, intended purpose

- Coverage of second- and third-degree skin loss due to burns or trauma
- Coverage of chronic non-healing wounds, such as venous, diabetic or pressure ulcers
- Coverage of autologous grafts
- Coverage of donor sites

Patient target group

Xenoderm is intended for use in patients with skin loss resulting from burns, trauma, chronic non-healing wounds, autologous grafts or donor sites, provided that no contraindications are present.

Use in pregnant or breastfeeding women, children and elderly patients is subject to an individual medical benefit-risk assessment by the treating physician.

Intended users

Xenoderm is intended to be used exclusively by medically trained healthcare professionals. Application must be performed under suitable aseptic conditions and in accordance with recognised medical and wound care practice.

Clinical benefit and performance characteristics

Xenoderm serves as a temporary biological skin replacement for covering and protecting wound areas. The product supports a moist wound environment, reduces fluid and heat loss, protects the underlying tissue against drying and external contamination, and may support granulation and re-epithelialisation as part of the natural wound healing process. Xenoderm has the following performance characteristics in particular:

- simple application after rehydration in sterile isotonic saline solution
- adapts closely to the wound surface with good adherence
- reduction of pain by covering exposed nerve endings
- no signs of vascularisation within the first seven days

- clinically non-antigenic effect
- protects the underlying tissue against bacterial ingress
- protects against drying
- dressing changes are usually carried out without pain or blood loss
- its transparency allows observation of the wound bed
- Xenoderm facilitates movement, early mobilisation and rehabilitation

Application

1. After sterile removal from the packaging, rehydration is performed at room temperature in sterile isotonic saline solution. The dressing should be soaked for at least 2 minutes or until the dressing becomes translucent.
2. At the time of initial or repeated application of Xenoderm, the wound area should be thoroughly debrided and gently cleansed, for example with an antiseptic. On second-degree wounds, Xenoderm may be left in place once adherent until it detaches during re-epithelialisation.
3. Press Xenoderm firmly onto the wound and smooth it out to avoid blister formation. Xenoderm can be cut to the required size using sterile scissors or another suitable sterile instrument and adapted to the wound margins.
4. As a precaution, Xenoderm should not be applied to weight-bearing areas of the body due to the risk of maceration.
5. Xenoderm can be fixed in place using gauze bandages or other suitable dressing materials.
6. Renewal of the Xenoderm dressing should be carried out according to the usual ward routine for dressing changes.

7. Xenoderm dressing changes are practically painless. Always moisten the dressing with sterile isotonic saline solution before removal.
8. The duration for which Xenoderm remains on the wound is variable and is determined by the course of wound healing.
9. Infected wounds with heavy exudation initially require more frequent dressing changes.
10. Daily Xenoderm dressing changes may be necessary for the above-mentioned wounds.
11. On all aseptic wounds, including granulating wounds, Xenoderm should be removed after no later than 10 days in order to prevent excessive adherence.

Contraindications

Xenoderm must not be used in the following cases:

- third-degree burns before removal of the eschar
- heavily infected wounds before appropriate cleansing or debridement
- patients with sensitisation or known hypersensitivity to materials of porcine origin

Interactions / combination with other preparations

Xenoderm must not be used in direct contact with ointments, powders or silicone preparations.

Xenoderm must not be brought into direct contact with methyl methacrylate.

Xenoderm should not be used in combination with antiseptics that release chlorine. Protein-damaging substances such as tannic acid or silver nitrate should also not be used.

Use during pregnancy and breastfeeding, use in children and elderly patients

No studies are available on the use of Xenoderm during pregnancy and breastfeeding or on its influence on human reproductive capacity.

Before using Xenoderm, the treating physician must therefore carry out an individual assessment of the benefit for the mother and the possible risks for the child.

There are no findings suggesting that special precautionary measures are necessary depending on the age of the patients to be treated.

Warnings, precautions and notes

- Xenoderm is sterilised by γ -irradiation and must not be resterilised.
- Xenoderm is individually packaged in sterile packaging.
- If the packaging is damaged, Xenoderm must be considered non-sterile and disposed of.
- Xenoderm has a shelf life of five (5) years when stored at room temperature.
- Xenoderm must not be used after the expiry date.
- Xenoderm must be protected from moisture.
- Xenoderm is a medical device and must be kept out of the reach of children.
- Xenoderm is intended for single use only.
- Xenoderm is intended exclusively for external, topical application as a temporary biological skin replacement on wound areas.
- Xenoderm must not be reused. After removal from the packaging and/or after contact with a patient, reuse carries an increased risk of contamination with subsequent risk of infection.

Disposal

Used or contaminated Xenoderm dressings must be handled and disposed of in accordance with recognised medical practice and the applicable national, regional and facility-specific regulations for medical or biological materials.

Non-contaminated packaging materials must be disposed of in accordance with applicable local regulations.

Where permitted by local regulations, used Xenoderm dressings may be disposed of with normal hospital waste.

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