

Hemostat Line



## Diacoll

Collagen wound dressing for fast and efficient hemostasis

**Properties:**

- Highly effective local hemostatic agent for wound care after venipuncture, e.g. dialysis.
- Shortens time to wound closure and minimizes risk of infection.
- Easy and rapid handling.
- Snuggles to the wound bed and can be removed easily (if required)

# Diacoll

CE 0499



Diacoll is an effective hemostatic wound dressing made from fully resorbable collagen.

The collagen matrix is triggering a fast, reliable and sustainable blood clotting. Due to easy handling Diacoll supports efficient hemostatic wound care, e.g. after dialysis. Diacoll minimizes treatment time for the patient and enables higher patient throughput.

### Further Properties:

- Class III medical product according to European guideline 93/42/EWG.
- Certified according to DIN EN ISO 13485:2012.
- Production following good manufacturing practice guidelines (GMP).
- Outstanding biocompatibility.
- Gamma sterilized.
- Fully resorbable.

### Specifications

Art-Nr.	Size	Contents
REF D0202K	20x20x7 mm	12 pieces per box
REF D3232K	32x32x12 mm	5 pieces per box



Health  
made in Germany