



Hemostat Line



# Surgicoll

Resorbable and highly effective hemostatic agent based on collagen

## Properties:

- Highly effective local hemostatic agent.
- Promotes wound healing.
- Fully resorbable within 1-2 weeks.
- Easy and simple handling.
- Available as fleece (Surgicoll), cone (Surgicoll Plug) or granulated powder (Surgicoll Cloud)

# Surgicoll

CE 0499



Surgicoll is a fully resorbable sponge like collagen matrix for effective hemostasis. Collagen is well known for its outstanding properties as hemostatic agent documented in numerous clinical reports.

The sponge like collagen matrix promotes reliable blood clotting due to two complementarily acting mechanisms. While the collagen matrix swelling properties cause a concentration of blood platelets, thus promoting blood clotting, the chemical and topological nature of collagen fibrils induces the release of blood clotting factors from aggregated platelets, thus triggering the coagulation cascade.

### Further Properties:

- Class III medical device according to European guideline 93/42/EWG.
- Certified according to DIN EN ISO 13485:2012.
- Production following good manufacturing practice guidelines (GMP).
- Outstanding biocompatibility.
- Gamma sterilized.
- Double sterile packaging for hygienic transfer to the surgery area.

### Specifications

| Product   | REF No.   | Size         | Content/Box |
|-----------|-----------|--------------|-------------|
| Surgicoll | S0101K    | 1x1x0,5 cm   | 10x3 pieces |
| Surgicoll | S0505K    | 5x5x0,5 cm   | 10 pieces   |
| Surgicoll | S0510K    | 5x10x0,5 cm  | 10 pieces   |
| Surgicoll | S1010K    | 10x10x0,5 cm | 5 pieces    |
| S. Cloud  | 1 025 PSC | 0,25 g       | 5 pouches   |
| S. Cloud  | 1 050 PSC | 0,50 g       | 5 pouches   |
| S. Plug   | ConeSP    | 0,01 g       | 5 pieces    |

