




Burns Line



# Xenoderm

Temporary biological skin replacement

## Properties:

- Protects the biological wound milieu.
- Promotes wound healing.
- Minimizes patient pain.
- Provides a microbial and mechanical barrier.
- Prevents electrolyte and heat loss.
- Excellent adherence to the wound bed.
- Moisturized removal with minimal pain.

# Xenoderm

CE 0499



Xenoderm is a split-skin derived acellular collagen matrix.

Large scale loss of skin leads to various complications and disturbs wound healing. The use of Xenoderm establishes a biological wound milieu promoting rapid epithelialization and good cosmetic results.

## Xenoderm

- Is designed as a semi-permanent cover for burns (degree II-III), traumatic and all full and partial thickness wounds.
- Is freeze dried dermis, harvested from animals destined for the human food chain.
- Is used where Homograft or Cadaver skin is used.
- Can be used on third degree burns immediately after excision where grafting is delayed for whatever reason.
- Can be the second component involved in the sandwich technique: **“Autograft and Cadaver Skin“** or **Xenoderm**.
- Can be applied on burns already covered with autograft, meshed or unmeshed and/or any skin substitute.
- Can be used in all trauma cases involving tissue loss or trauma cases where wound closure is delayed.
- Can be used to cover non healing wounds as a temporarily wound dressing, alone or in combination with others.
- Is close to human skin in its biocompatibility and bears no risk of human virus transmission.

## Further Properties:

- Class III medical product according to European Guideline 93/42/EWG.
- Manufactured in accordance with a quality management system certified to DIN EN ISO 13485:2012.
- Production following good manufacturing practice guidelines (GMP).
- Gamma sterilized.

## Specifications

REF-No.	Size	Contents
REF 0505X	50x50x0,3 mm	10 pieces per box
REF 1010X	100x100x0,3 mm	5 pieces per box
REF 1020X	100x200x0,3 mm	5 pieces per box
REF 1030X	100x300x0,3 mm	5 pieces per box
REF 1040X	100x400x0,3 mm	5 pieces per box

